NTI-Kahla GmbH

Rotary Dental Instruments

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When using hand type instruments please put on the safety chain for root canal instruments or work with rubber dam. When using a hand type endo instrument in an endo contra-angle (Giomatik contra-angle), please follow the manufacturer's instructions.

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K reamer, No. 171... :

Exert pressure during use and rotate clockwise up to 90°.

For power-assisted canal preparation, use a contra-angle with yellow colour-coding -1/4 turn. Speed: 450 - 800 r.p.m.

K files, No. 173... :

Filing motion combined with clockwise rotary movements up to 45°.

For power-assisted canal preparation, use a contra-angle with yellow colour-coding - 1/4 turn. Speed: 450 - 800 r.p.m.

Hedstroem files, No. 174... :

To avoid distorting handtype instruments, they should only be used with reciprocating movements. For power-assisted canal preparation, use a contra-angle with yellow colour-coding - up to 1/4 turn.

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To prevent the instrument drying out after use, immerse it in a tank full of disinfectant approved by the German Society for Hygiene and Microbiology (DGHM).

(Always follow the manufacturer's recommended mixing ratio and immersion time in the disinfectant. Do not use disinfectants containing phenic acid). Rinse the instrument under running water and check for residues, clean it again if necessary and dry it.

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Phone: +49(36424)573-0
Fax: +49(36424)573-29
http://www.nti.de

GB Instructions for endodontic instruments:

When using hand type instruments please put on the safety chain for root canal instruments or work with rubber dam. When using a hand type endo instrument in an endo contra-angle (Giromatik contra-angle), please follow the manufacturer's instructions.

Peeso, No. 183.15... and 183.19... :

Available in two lengths: 15 mm and 19 mm. Use with a contra-angle with green colour-coding at 800 to 1,200 r.p.m.

Gates, No. 180.15... and 180.19... :

Available in two lengths: 15 mm and 19 mm. Use with a contra-angle with green colour-coding at 450 to 800 r.p.m.

Reamer, No. 183.20... :

Use with a contra-angle with green colour-coding at 800 to 1,200 r.p.m.

Barbed broaches, No. 9107... :

Insert the broach into the root canal and rotate it through 180° - 360° to engage and withdraw the pulp.

K reamer, No. 171... :

Exert pressure during use and rotate clockwise up to 90°. For power-assisted canal preparation, use a contra-angle with yellow colour-coding -1/4 turn. Speed: 450 - 800 r.p.m.

K files, No. 173... :

Filing motion combined with clockwise rotary movements up to 45°. For power-assisted canal preparation, use a contra-angle with yellow colour-coding - 1/4 turn. Speed: 450 - 800 r.p.m.

Hedstrom files, No. 174... :

To avoid distorting handtype instruments, they should only be used with reciprocating movements. For power-assisted canal preparation, use a contra-angle with yellow colour-coding - up to 1/4 turn. Speed: 450 - 800 r.p.m.

Root fillers, Nos. 178... and 178S... :

Use a green colour-coded contra-angle. Immerse the spiral (Lentulo) in the root canal sealer and, with the handpiece at a standstill, insert it as deep as possible into the root canal. Rotate the spiral at low speed (approx. 800 r.p.m.) and slowly withdraw it from the canal. The sealer remains in the canal. Product No. 178S... , Safety spiral filler on a contra-angle shank - can be inserted into the root canal flexibly.

Using Instruments for contra-angles:

When using a contra-angle (power-assisted preparation), always take care and follow the manufacturer's instructions. Incorrect operation can easily cause the instrument to fracture.

Instrument coding:

Hand types and contra-angle shanks are colour-coded to denote their size and standard. Hand type instruments are also marked with a symbol and size. The plastic grips of the barbed broaches are coloured to denote their size. Spiral fillers with metal shanks are coded with coloured rings.

Servicing and storing the instruments hygienically:

To prevent the instrument drying out after use, immerse it in a tank full of disinfectant approved by the German Society for Hygiene and Microbiology (DGHM). (Always follow the manufacturer's recommended mixing ratio and immersion time in the disinfectant. Do not use disinfectants containing phenic acid). Rinse the instrument under running water and check for residues, clean it again if necessary and dry it.

Sterilization:

	Temperature	Pressure	Time
Steam (autoclave)	134°C	2,1 bar	5 mins.
	121°C	1,1 bar	20 mins.
Hot-air sterilizer	180°C		30 mins.

Storage and Keeping Conditions:

In dry conditions and protected against contaminants.

Paperpoints: For drying root canal. Sterility is not guaranteed after the pack is opened.

Guttaperchapoints: For root canal filling. Do not autoclave or heat sterilize. Desinfect before use.

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