NTI-Kahla GmbH

Rotary Dental Instruments

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TRI-STAR, TRI-JET

DESCRIPTION

Parapulpal retention pins, for use by dentist only

TYPES

See trade names

COMPOSITION

* Implant alloy => titanium alloy

* Pins: Ti 6Al 4V

* According to the literature, material considered hypoallergenic and biocompatible

* (Goldplated) steel pins = stainless steel 1.4305

INDICATIONS

Retention pins for the build-up of fractured teeth, corners; twist anchorage

CONTRAINDICATIONS

Insufficient healthy, caries-free dentine

SIDE EFFECTS

None known, except for metal allergy

PRECAUTIONS

1. Avoid contact with pulp.
2. Sterilise BEFORE use.
3. Observe general rules of hygiene.
4. Preserve distance of approx. 1 mm to dental enamel.
5. Use appropriate bender.
6. In case of composite caps, use veneer.

INTERACTIONS None known

APPLICATION

- ① Pre-drill until the depth-limit
- ② Insert pin using contra angle handpiece or insert pin using blue handle
- ③ Pins in situ
- ④ Anatomical pin adjustment
- ⑤ Build-up and restoration completed

STORAGE INFORMATION: Dry

Mandrels are reusable! After inserting the pin remove the remainder from the mandrel by twisting it out. Then simply screw in a new pin from the refill-box. Cleaning: Mandrels are made of anodized aluminium. Use only appropriate cleaning solutions for aluminium, or clean under running water and finally sterilise in autoclave.

PRESENTATIONS:

Standard Assortment of 25 pins, 2 twist drill, 1 handle;

Economy Assortment: 50 pins, 2 twist drills, 1 handle;

Refill pack of 100 pins.

HYGIENE INFORMATION: Autoclave: 134 °C - 2.1 bar- 5 minutes

LIABILITY: Within the framework of the customary material defects guarantee, no liability whatsoever can be accepted for improper use.

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