
NTI-Kahla GmbH • Rotary Dental Instruments
Im Camisch 3 • D-07768 Kahla/Germany
Tel. +49(0)36 424-573-0 • Fax +49(0)36 424-573-29
e-mail: nti@nti.de • http://www.nti.de

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GB Product Information FiberMaster Glass fibre root posts

Customer information

Please read the product information carefully before using FiberMaster posts. The FiberMaster material was specially developed for use in dentistry and should be used in accordance with the recommended technique. The manufacturer does not accept any liability for damage caused by disregarding these instructions or incorrect use. The user is personally responsible for checking that the material is suitable for the intended purpose, especially if this is not included in the recommended uses.

Contents:

1. Application of the system

1.1. Requirements for using FiberMaster posts

2. Recommended technique

3. Disinfection and cleaning

1. Application of FiberMaster posts

Teeth with badly damaged crowns or fractured teeth with no other lesions of the roots and still completely intact periodontally

1.1. Requirements for using the FiberMaster system

- The root should be periodontally intact.
- The lesion should not extend below the gingival margin.
- The residual tooth structure should be caries-free.
- Successfully completed root treatment.
- The post should fit against the root walls.
- Anchorage depth at least 2/3 of the length of the root canal.
- Adequate root wall thickness (min. 1.5 mm).
- As large a post diameter as possible.
- The occlusal or incisal surface of the prepared tooth should be prepared for a core.
- Preparation margin at least 1 mm into the natural tooth structure.

2. Recommended technique

- Select the size of the root post using the x-ray.
- Decide at this stage whether to use a post with or without a head. Several posts with or without a head can be used for teeth with multiple roots, depending on the space available.
- Remove the root filling.
- Use a root-facing bur to prepare a rest for the FiberMaster head by cutting to a depth of approx. 1 mm (2,000 rpm) on the surface of the prepared tooth. This stage is not required for FiberMaster posts without a head.
- Prepare the root canal to the required depth with a drill. The correct depth is determined on the x-ray. Stoppers from the endodontic range can be put on the drill to mark the required length (motor speed: 500 - 1,000 rpm). When drilling the root canal, ensure that a minimum of approx. 4 mm root filling is left in the apex.
- FiberMaster posts can be easily cut to the appropriate length using a diamond disc (only use discs without holes or slits). The post should be shortened extraorally! The posts are radiopaque, so the fit can be checked on the x-ray.
- Disinfect the post with 70% medical alcohol.

The posts cannot be sterilised as the temperatures damage the material.

- The alcohol is also used to remove any grease (hands) to ensure that the adhesive bonds optimally. Do not handle after cleaning.
- Use an adhesive luting technique. A rubber dam should always be used with this technique. **Adhesive technique:**
 - etch the root canal and surface of the prepared

tooth with acid

- prepare the dentine walls and poss.

FiberMaster post with primer

- apply bonder if required.

Apply luting composite to the root post.

- Use a low viscosity dual composite for adhesive. High viscosity dual composites should only be used with an ultrasonic unit. When cementing the post, apply light pressure to the occlusal of the post with a ball instrument. Dual composites should be used as they harden without light curing. Adhere strictly to the instructions for use of the adhesive manufacturer! Use a rubber dam and ensure the post cannot be swallowed or aspirated.
- **Please note:**
Do not clean the root canal with NaOCl (sodium hypochloride) or H₂O₂ (hydrogen peroxide), as these chemicals can have a detrimental effect on the polymerisation of the composite. Use paper tips for drying. Compressed air may contain oil which weakens the bond.
- Insert the root post into the canal applying gentle pressure to allow any excess to escape. Carefully remove the excess and light cure. Caution: The root post cannot transmit light to the required depth. This is why self-curing or dual composite adhesive should be used.
- Fabricate the core using a suitable core composite. (Adhere to the manufacturer's instructions).
- Complete the preparation ensuring that it is at least 1 mm below the root surface.
- If necessary, the size or shape of the head can easily be adjusted intraorally using a diamond instrument. A rubber dam should be used when adjusting.
- The posts should be secured with a safety chain or other safety measures taken to prevent aspiration or swallowing by the patient.

3. Désinfection and cleaning

The rotary instruments are disinfected and sterilised in the usual manner. The posts are disinfected with 70% medical alcohol and any grease cleaned off.

The posts cannot be sterilised.

Do not handle the posts after disinfection.

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© 11/2006, NTI-Kahla GmbH Rotary Dental Instruments